

ORTHOPAEDIC

Bone cutters
Bone Rongeurs
Bone Curettes
Bone Holding Forceps
Bone Holding Clamps
Plaster Spreaders and Shears
Bone Levers
Retractors
Dissecting Forceps and Hooks
Raspatories, Burr and Rugines
Osteotomes, Chisels and Gouges
Wire Cutters and Pliers
Accessories
Mallets
Raspatories, Files and Punches
Pushers, Impactors and Skids
Elevators
Retractors and Bone Levers
Saws, Shears and Knives
Screw and Plate Instruments
Pins
Sterile Kirschner and Arthrodesis Wires
Extremities Instruments and Sets
Shoulder Instruments



+92-340-6112989




www.rizacosurgical.com

Bone Cutters

McIndoe Forceps

Lightweight, Double action, angled jaw with prismatic points



Code	Sizes
RS-OR/001	165mm long
RS-OR/002	191mm long

Ruskin Liston

Double action, straight, angled



Straight

Code	Sizes
RS-OR/003	140mm long
RS-OR/004	145mm long
RS-OR/005	190mm long

Angled

Code	Sizes
RS-OR/006	190mm long

Liston Screw Joint Forceps

Straight, angled



Straight

Code	Sizes
RS-OR/007	150mm long
RS-OR/008	180mm long
RS-OR/009	200mm long
RS-OR/010	230mm long

Angled

Code	Sizes
RS-OR/011	150mm long
RS-OR/012	180mm long

Liston Box Joint Forceps



Straight

Code	Sizes
RS-OR/013	140mm long
RS-OR/014	200mm long
RS-OR/015	220mm long

Angled

Code	Sizes
RS-OR/016	170mm long
RS-OR/017	220mm long

Liston double Action Forceps


Straight, for heavy duty




Code	Sizes
RS-OR/018	305mm long
RS-OR/019	300mm long
RS-OR/020	190mm long

Boehler Bone Forceps

Double action, Straight, angled




Code	Sizes
RS-OR/021	145mm long
RS-OR/022	150mm long
RS-OR/023	170mm long
RS-OR/024	190mm long



Cleveland Bone Forceps


Curved, angled




Code	Sizes
RS-OR/025	150mm long
RS-OR/026	170mm long

Stille-Horsley Forceps

Double action, Straight, angled




Code	Sizes
RS-OR/027	255mm long




Stille-Liston 90 degree jaw Forceps

Double action, 90 degree jaw





Code	Sizes
RS-OR/028	267mm long



Stille-Liston Forceps



Double action, angled on flat, straight



Code	Sizes
RS-OR/029	180mm long
RS-OR/030	270mm long

Ruskin-Roland Forceps

Double action, angled on flat, straight




Code	Sizes
RS-OR/031	171mm long
RS-OR/032	184mm long

Bone Rongeurs

Adson Rongeur

Double action, Straight, angled




Code	Sizes
RS-OR/033	200mm long x 7mm tip
RS-OR/034	200mm long x 8mm tip

7mm 8mm

Bacon Cranial Rongeur


Curved, angled



Code	Sizes
RS-OR/035	200mm long x 5mm tip

Beyer Rongeur

Double action, Straight, angled




Code	Sizes
RS-OR/036	178mm long x 3mm tip
RS-OR/037	178mm long x 4mm tip

3mm 4mm

Blumenthal Rongeur

Double action, 90 degree jaw




Code	Sizes
RS-OR/038	152mm long - 30° angle
RS-OR/039	152mm long - 45° angle
RS-OR/040	152mm long - 90° angle

30° 45° 90°

Cleveland Rongeur

Curved




Code	Sizes
RS-OR/041	140mm long
RS-OR/042	170mm long

3mm 4mm

Dale Rongeur

Double action, curved




Code	Sizes
RS-OR/043	368mm long
RS-OR/044	229mm long

3mm
4mm
5mm

DeVilbiss Cranial Rongeur


Single action



Code	Sizes
RS-OR/045	210mm long x 7mm tip

Echlin Duckbill Rongeur


Double action, angular



Code	Sizes
RS-OR/051	230mm long x 2mm tip
RS-OR/052	230mm long x 3mm tip
RS-OR/053	230mm long x 4mm tip

Hartmann Rongeur


Single action, sharp, curved jaws



Code	Sizes
RS-OR/057	185mm long x 4mm tip

Friedman Rongeur

Single action




Straight

Code	Sizes
RS-OR/046	140mm long x 2mm tip
RS-OR/047	140mm long x 3mm tip

Angled


Code	Sizes
RS-OR/048	140mm long x 3mm - slight curve
RS-OR/049	140mm long x 3mm - full curve
RS-OR/050	140mm long x 4mm - slight curve

Straight



Friedman Rongeur Micro

Straight, Slight curved, curved 90



Code	Sizes
RS-OR/054	140mm long x 1.7mm tip (a)
RS-OR/055	140mm long x 1.7mm tip (b)
RS-OR/056	140mm long x 1.7mm tip (c)

custom sizes

Jansen Rongeur

single action, straight, curved




Straight

Code	Sizes
RS-OR/058	178mm long
RS-OR/059	190mm long


Curved

Code	Sizes
RS-OR/060	178mm long
RS-OR/061	190mm long



Juers-Lempert Rongeur

Single action, curved handle




Code	Sizes
RS-OR/062	184mm long x 3.5mm tip

Leksell Rongeur

Double action, full curved, 240mm length




Code	Sizes
RS-OR/067	3mm jaw
RS-OR/068	4mm jaw
RS-OR/069	5mm jaw
RS-OR/070	6mm jaw
RS-OR/071	7mm jaw
RS-OR/072	8mm jaw



Leksell-Stille Rongeur


Double action, curved, 240mm length



Code	Sizes
RS-OR/074	3mm jaw
RS-OR/075	4mm jaw
RS-OR/076	5mm jaw
RS-OR/077	6mm jaw
RS-OR/078	7mm jaw
RS-OR/079	8mm jaw

Kleinert-Kutz Rongeur

Double action, full or slight curved




Straight

Code	Sizes
RS-OR/063	152mm long x 2mm tip
RS-OR/064	152mm long x 3mm tip


Angled

Code	Sizes
RS-OR/065	152mm long x 2mm tip
RS-OR/066	152mm long x 3mm tip



Leksell-Stille Rongeur

Double action, slightly angled jaw



Code	Sizes
RS-OR/073	240mm L x 8mm jaw

Pair Replacement Springs for Leksell-Stille Rongeur


RS-OR/067 to RS-OR/079



Code	Sizes
RS-OR/080	For Leksell-stille 240mm

Lempert Rongeur

Single action action, Curved 3mm




Straight

Code	Sizes
RS-OR/081	152mm long



Angled

Code	Sizes
RS-OR/082	152mm long



Lempert Rongeur

Single action action, Curved 3mm



Straight

Code	Sizes
RS-OR/083	190mm long




Angled

Code	Sizes
RS-OR/084	190mm long



Marquardt Rongeur

Double action, Slight curved 3mm



Code	Sizes
RS-OR/085	200mm long



Mead Rongeur

Single action, curved 4mm




Code	Sizes
RS-OR/086	165mm long



Mayfield Rongeur

Double action, 3mm Tip




Code	Sizes
RS-OR/087	175mm long




Ortho Rongeur

Ferris Smith Kerrison handle, 254mm long




Code	Sizes
RS-OR/088	10 x 16mm
RS-OR/089	7 x 16mm
RS-OR/090	5 x 14mm

Tip Sizes



Rat Tooth Rongeur

single action, Tip 8mm




Code	Sizes
RS-OR/091	200mm long



Reiner Rongeur

Single action action, 178mm long



Tip Sizes

Code	Sizes
RS-OR/092	3mm curved
RS-OR/093	4mm curved

Ruskin Rongeur Curved

Double action, 184mm length



Tip Sizes

Code	Sizes
RS-OR/094	4mm
RS-OR/095	5mm
RS-OR/096	6mm
RS-OR/097	7mm




Ruskin Rongeur Straight

double action, 184mm length





Tip Sizes

Code	Sizes
RS-OR/098	4mm
RS-OR/099	5mm
RS-OR/100	6mm
RS-OR/101	7mm



Ruskin Rongeur 5mm Curved


Double action, 5mm Tip



Code	Sizes
RS-OR/102	184mm long

Sauerbruch Rongeur

Double action, 300mm long



Tip Sizes

Code	Sizes
RS-OR/103	300mm Square
RS-OR/104	300mm Round

Semb Rongeur

Double action, 17mm Round tip



Code	Sizes
RS-OR/105	225mm long



Smith Peterson Rongeur

Double action, 241mm long



Straight


Code	Sizes
RS-OR/106	3mm curved

Curve

Code	Sizes
RS-OR/107	3mm curved

Full curve

Code	Sizes
RS-OR/108	3mm curved




Stellbrink Rongeur

Single action, 165 to 175mm length
Lap Joint



Tip Sizes

Code	Sizes
RS-OR/109	2mm



Stellbrink Rongeur

Single action, 171mm length



Slight curve

Code	Sizes
RS-OR/110	2mm


Full curve

Code	Sizes
RS-OR/111	2mm



Stille Ruskin Rongeur

Double action, Curved 8mm Tip



Code	Sizes
RS-OR/112	240mm long

Stille-Luer Rongeur

Double action, straight, curved



26.5cm length x 17mm tip

Code	Sizes
RS-OR/113	Straight
RS-OR/114	Curved

26.5cm length x 10mm tip

Code	Sizes
RS-OR/115	Straight
RS-OR/116	Curved

24.1cm length x 17mm tip

Code	Sizes
RS-OR/117	Straight
RS-OR/118	Curved



Stille-Luer Rongeur

Double action



21.5cm length x 10mm tip

Code	Sizes
RS-OR/119	Straight
RS-OR/120	Curved

Stille-Luer Rongeur Angular


Double action, 228mm long



Code	Sizes
RS-OR/121	5mm
RS-OR/122	6mm

Stookey Cranial Rongeur

Single action, Angled




191mm Length

Code	Sizes
RS-OR/123	3mm

Sypert Rongeur

Double action, 368mm length = 14½"




Tip sizes

Code	Sizes
RS-OR/124	3mm
RS-OR/125	4mm
RS-OR/126	5mm
RS-OR/127	8mm
RS-OR/128	12mm

Sypert Rongeur with Teeth

Double action, Curved 8mm Tip
368mm length = 14½"



Tip sizes

Code	Sizes
RS-OR/129	5mm
RS-OR/130	8mm
RS-OR/131	12mm

Zaufal-Jansen Rongeur

Double action, curved, 180mm length



Tip sizes

Code	Sizes
RS-OR/132	4mm

Orthopedic Currettes Bone

Bone Curette

Hollow Handle, straight,
381mm long, oval Tips

15 inches long



Tip sizes

Code	Sizes
RS-OR/133	(6.9mm) #1
RS-OR/134	(6.0mm) #1/0
RS-OR/135	(7.8mm) #2
RS-OR/136	(5.1mm) #2/0
RS-OR/137	(8.7mm) #3
RS-OR/138	(4.2mm) #3/0
RS-OR/139	(9.6mm) #4
RS-OR/140	(3.3mm) #4/0
RS-OR/141	(10.5mm) #5
RS-OR/142	(11.4mm) #6



Tip



straight

Bone Curette

Hollow Handle, 30 degree angled,
381mm long, oval Tips

15 inches long



Tip sizes

Code	Sizes
RS-OR/143	(6.9mm) #1
RS-OR/144	(6.0mm) #1/0
RS-OR/145	(7.8mm) #2
RS-OR/146	(5.1mm) #2/0
RS-OR/147	(8.7mm) #3
RS-OR/148	(4.2mm) #3/0
RS-OR/149	(9.6mm) #4
RS-OR/150	(3.3mm) #4/0
RS-OR/151	(10.5mm) #5
RS-OR/152	(11.4mm) #6



Tip



30 degree
angled

Bone Curette

straight, double handed knurled T-handle,
381mm long

15 inches long



Tip sizes

Code	Sizes
RS-OR/153	11x20mm
RS-OR/154	15x28mm
RS-OR/155	18x34mm
RS-OR/156	8.4x18.2mm



Brun Curette 6 1/4"


Hollow handle straight, 158mm 1/4 length,
Round

6 1/4 inches long



Tip sizes

Code	Sizes
RS-OR/157	(6.9mm) #0
RS-OR/158	(6.0mm) #1
RS-OR/159	(7.8mm) #2
RS-OR/160	(5.1mm) #2/0
RS-OR/161	(8.7mm) #3
RS-OR/162	(4.2mm) #3/0
RS-OR/163	(9.6mm) #4
RS-OR/164	(3.3mm) #4/0
RS-OR/165	(10.5mm) #5
RS-OR/166	(11.4mm) #6



Tip

Brun Curette 7"

Hollow handle, angled, oval
178mm length



Tip sizes

Code	Sizes
RS-OR/167	#0
RS-OR/168	#1
RS-OR/169	#2
RS-OR/170	#2/0
RS-OR/171	#3
RS-OR/172	#3/0
RS-OR/173	#4
RS-OR/174	#4/0
RS-OR/175	#5
RS-OR/176	#5/0
RS-OR/177	#6



Angled

Brun Curette 7"

Hollow handle, straight, oval,
178mm length



Tip sizes

Code	Sizes
RS-OR/178	#0
RS-OR/179	#1
RS-OR/180	#2
RS-OR/181	#2/0
RS-OR/182	#3
RS-OR/183	#3/0
RS-OR/184	#4
RS-OR/185	#4/0
RS-OR/186	#5
RS-OR/187	#5/0
RS-OR/188	#6



Straight

Brun Curette 8"

Hollow handle, angled, oval
203mm length



Tip sizes

Code	Sizes
RS-OR/189	#0
RS-OR/190	#1
RS-OR/191	#2
RS-OR/192	#2/0
RS-OR/193	#3
RS-OR/194	#3/0
RS-OR/195	#4
RS-OR/196	#4/0
RS-OR/197	#5
RS-OR/198	#5/0
RS-OR/199	#6



Brun Curette 8"

Hollow handle, straight, oval
203mm length



Tip sizes

Code	Sizes
RS-OR/200	#0
RS-OR/201	#1
RS-OR/202	#2
RS-OR/203	#2/0
RS-OR/204	#3
RS-OR/205	#3/0
RS-OR/206	#4
RS-OR/207	#4/0
RS-OR/208	#5
RS-OR/209	#5/0
RS-OR/210	#6



Brun Curette 9"

Hollow handle, angled, oval
228mm length



Tip sizes

Code	Sizes
RS-OR/211	#0
RS-OR/212	#1
RS-OR/213	#2
RS-OR/214	#2/0
RS-OR/215	#3
RS-OR/216	#3/0
RS-OR/217	#4
RS-OR/218	#4/0
RS-OR/219	#5
RS-OR/220	#5/0
RS-OR/221	#6



Brun Curette 9"

Hollow handle, straight, oval
228mm length



Tip sizes

Code	Sizes
RS-OR/222	#0
RS-OR/223	#1
RS-OR/224	#2
RS-OR/225	#2/0
RS-OR/226	#3
RS-OR/227	#3/0
RS-OR/228	#4
RS-OR/229	#4/0
RS-OR/230	#5
RS-OR/231	#5/0
RS-OR/232	#6



Brun Curette 9"

Hexagonal handle, angled, oval



Tip sizes

Code	Sizes
RS-OR/233	#0
RS-OR/234	#1
RS-OR/235	#2
RS-OR/236	#2/0
RS-OR/237	#3
RS-OR/238	#3/0
RS-OR/239	#4
RS-OR/240	#4/0
RS-OR/241	#5
RS-OR/242	#6



Brun Curette 9"

Hexagonal handle, Straight, oval



Tip sizes

Code	Sizes
RS-OR/243	#0
RS-OR/244	#1
RS-OR/245	#2
RS-OR/246	#2/0
RS-OR/247	#3
RS-OR/248	#3/0
RS-OR/249	#4
RS-OR/250	#4/0
RS-OR/251	#5
RS-OR/252	#6




Bushe Curette 10"

Hollow handle, angled, oval
254mm length



Tip sizes

Code	Sizes
RS-OR/253	(5.2mm)#0
RS-OR/254	(6.0mm)#1
RS-OR/255	(4.4mm)#2/0
RS-OR/256	(3.6mm)#3/0
RS-OR/257	3mm reverse bend
RS-OR/258	4mm reverse bend
RS-OR/259	5mm reverse bend



Bushe Curette 10"

Hollow handle, straight, oval
254mm length



Tip sizes

Code	Sizes
RS-OR/260	(5.2mm)#0
RS-OR/261	(6.0mm)#1
RS-OR/262	(4.4mm)#2/0
RS-OR/263	(3.6mm)#3/0




Caspar Bone Curette

Angled, toothed, sharp/blunt,
266mm length



Tip sizes

Code	Sizes
RS-OR/264	3mm
RS-OR/265	4mm




Caspar Bone Curette

Straight, toothed, sharp/blunt,
254mm length



Tip sizes

Code	Sizes
RS-OR/266	4mm
RS-OR/267	5mm



Cobb Curette

Angled, 279mm Length, aluminum handle,
stainless



Tip sizes

Code	Sizes
RS-OR/268	#0
RS-OR/269	#1
RS-OR/270	#2
RS-OR/271	#2/0
RS-OR/272	#3
RS-OR/273	#3/0
RS-OR/274	#4
RS-OR/275	#4/0
RS-OR/276	#5
RS-OR/277	#6




Cobb Curette

Straight, 279mm Length, aluminum handle,
stainless



Tip sizes

Code	Sizes
RS-OR/278	#0
RS-OR/279	#1
RS-OR/280	#2
RS-OR/281	#2/0
RS-OR/282	#3
RS-OR/284	#3/0
RS-OR/285	#4
RS-OR/286	#4/0
RS-OR/287	#5
RS-OR/288	#6



Cone Ring Curette

25 degree angled, aluminum handle, sharp/sharp, 381mm length



Tip sizes

Code	Sizes
RS-OR/289	(3mm)#1
RS-OR/290	(6mm)#2
RS-OR/291	(8mm)#3

Cone Ring Curette

Straight, aluminum handle, sharp/sharp, 381mm length



Tip sizes

Code	Sizes
RS-OR/292	(3mm)#1
RS-OR/293	(6mm)#2
RS-OR/294	(8mm)#3

Cone Ring Curette 9"

25 degree angled, aluminum handle, sharp/sharp



Tip sizes

Code	Sizes
RS-OR/295	(3mm)#1
RS-OR/296	(6mm)#2
RS-OR/297	(8mm)#3

Cone Ring Curette 9"

Straight, aluminum handle, sharp/sharp



Tip sizes

Code	Sizes
RS-OR/298	(3mm)#1
RS-OR/299	(6mm)#2
RS-OR/300	(8mm)#3


Cone Ring Curette 17"

Angled, 430mm Length, phenolic handle sharp/sharp



Tip sizes

Code	Sizes
RS-OR/301	3mm




Cone Ring Curette 20"

Angled, 508mm Length, phenolic handle sharp/sharp



Tip sizes

Code	Sizes
RS-OR/302	6mm



Cup Curette

Forward angle, black handle,
12 ½" 13 ⅓" length



12 ½" Tip sizes

Code	Sizes
RS-OR/303	7mm
RS-OR/304	9mm
RS-OR/305	11mm

13 ⅓" Tip sizes

Code	Sizes
RS-OR/306	7mm
RS-OR/307	9mm
RS-OR/308	11mm



Cup Curette

Reverse angle, black handle,
12 ½" 13 ⅓" length



12 ½" Tip sizes

Code	Sizes
RS-OR/321	7mm
RS-OR/322	9mm

13 ⅓" Tip sizes

Code	Sizes
RS-OR/323	7mm
RS-OR/324	9mm



Halle Curette

Straight, 215mm Length- 8½" inches,
hollow handle, oval tip



8 ½" Tip sizes

Code	Sizes
RS-OR/335	4.8mm
RS-OR/336	5.8mm
RS-OR/337	6.8mm



Cup Curette

Straight, black handle,
12 ½" 13 ⅓" length



12 ½" Tip sizes

Code	Sizes
RS-OR/309	7mm
RS-OR/310	11mm
RS-OR/311	13mm
RS-OR/312	15mm
RS-OR/313	17mm
RS-OR/314	19mm

13 ⅓" Tip sizes

Code	Sizes
RS-OR/315	7mm
RS-OR/316	11mm
RS-OR/317	13mm
RS-OR/318	15mm
RS-OR/319	17mm
RS-OR/320	19mm



Epstein Curette

Reverse angle, hollow handle,
203mm length, 8" inches



8" Tip sizes

Code	Sizes
RS-OR/325	#0
RS-OR/326	#1
RS-OR/327	#2
RS-OR/328	#2/0
RS-OR/329	#3
RS-OR/330	#3/0
RS-OR/331	#4
RS-OR/332	#4/0
RS-OR/333	#5
RS-OR/334	#6



Hardy Bayonet Curette

9 ½", 4 ¾ length, 45 degree right,
241mm x 120mm



9 ½" Tip sizes

Code	Sizes
RS-OR/338	3mm
RS-OR/339	5mm